



GE HealthCare

Versana Essential™

Ultrasound System Specification sheet

Product description

Versana Essential ultrasound helps you provide world-class care with no compromise on quality. This practical, versatile system has broad applicability and is well-suited for clinicians of all experience levels who wish to upgrade their capabilities. Versana Essential is a good fit for private clinics, GP offices, OB/GYN offices and other primary care settings as it is affordable, enables comprehensive scanning capability with excellent image quality, and is built to perform reliably in busy care environments, day after day.



Contents

1. General specifications	3	
1.01 Dimensions and weight	3	
1.02 Electrical power	3	
1.03 Control design	3	
2. User interface	3	
2.01 Operator keyboard	3	
2.02 Monitor	3	
3. System overview	3	
3.01 System specifications	3	
3.02 Applications	3	
3.03 Scanning methods	3	
3.04 Transducer type	4	
3.05 Operation modes	4	
4. System standard features	4	
5. System options	4	
6. Peripheral options	5	
7. Display modes	5	
7.01 Live and stored display format	5	
8. Selectable alternating modes	5	
9. Display annotation	5	
9.01 General user interface	5	
9.02 B-Mode	6	
9.03 Color Flow Mode	6	
9.04 PDI Mode	6	
9.05 M-Mode	6	
9.06 Doppler Mode	6	
9.07 Easy 3D	6	
9.08 Advanced 3D	6	
10. General system parameters	7	
10.01 System setup	7	
10.02 System scanning parameters	7	
10.03 B-Mode	7	
10.04 Coded Harmonic Imaging	8	
10.05 SRI-HD	8	
10.06 CrossXBeam	8	
10.07 Color Flow Mode	8	
10.08 Power Doppler Imaging Mode	9	
10.09 M-Mode	9	
10.10 Anatomical M-Mode (option)	9	
10.11 Pulse Wave Doppler Mode	9	
10.12 Continues Wave Doppler Mode	10	
10.13 Cine memory/image memory	10	
10.14 Image storage	10	
10.15 Connectivity and DICOM®	10	
10.16 Virtual convex	10	
10.17 LOGIQView (option)	11	
10.18 Easy 3D (option)	11	
10.19 Advanced 3D (option)	11	
10.20 TVI (option)	11	
10.21 TVM (option)	11	
10.22 Follow-up tool (option)	11	
10.23 Needle recognition (option)	11	
10.24 Scan Coach (option)	11	
10.25 My Trainer	11	
10.26 Battery (option)	11	
10.27 Scan Assistant	11	
10.28 RSVP	12	
10.29 Whizz	12	
10.30 Whizz CF mode	12	
10.31 Lateral gain compensation (LGC)	12	
10.32 Elastography (option)	12	
10.33 LI-RADS®ACR (Option)	12	
10.34 e-Delivery	12	
10.35 Probe Check	12	
10.36 Controls available while “live”	12	
10.37 Controls available on Freeze or Recall	12	
11. Measurements/calculations	13	
11.01 General B-Mode	13	
11.02 General M-Mode	13	
11.03 General Doppler measurements/calculations	13	
11.04 Real-time Doppler auto measurements/calculations	13	
11.05 OB measurements/calculations	14	
11.06 GYN measurements/calculations	14	
11.07 Vascular measurements/calculations	14	
11.08 Urological calculations	14	
12. Cardiac measurements/calculations	15	
12.01 B-Mode measurements	15	
12.02 M-Mode measurements	15	
12.03 Doppler Mode measurements	15	
12.04 Color Flow Mode measurements	17	
12.05 Combination Mode measurements	17	
12.06 Cardiac workstation	17	
13. Probes	18	
13.01 4C-RS	18	
13.02 L6-12-RS	18	
13.03 LK760-RS	18	
13.04 E8C-RS	18	
13.05 8C-RS	19	
13.06 3Sc-RS	19	
13.07 6S-RS	19	
13.08 E8Cs-RS	19	
14. Inputs and outputs	20	
15. Safety conformance	20	

1. General specifications

1.01 Dimensions and weight

Height with monitor	• 1395 mm (54.9 in)
Width	• Keyboard: 560 mm (22.0 in) • Caster: 520 mm (20.5 in)
Depth	• Maximum: 620 mm (24.4 in) • Caster: 610 mm (24.0 in)
Weight (no Peripherals)	Less than 45 kg

2.02 Monitor

Fixed monitor arm
Tilt/Rotate/Pan
Tilt angle: +25°/-90° Rotate angle: -90°, +90° Brightness/contrast/ color temperature adjustment Fold-down for transportation

1.02 Electrical power

Voltage 100 – 240 VAC

Frequency 50/60 Hz

Power consumption maximum of 350 VA with peripherals

1.03 Control Design

Max 3 active probe ports

Integrated SSD (256 GB)

Integrated speakers

Probe holders, removable for cleaning and washing

Gel holder, removable for cleaning and washing

Front and rear handles

Probe cable management slots

Easily removable air filters

Wheels: Locking mechanism that provides rolling lock and
caster swivel lock

Low Noise Design: less than 28dB (measured at 50cm far
from console, 25°C room temperature)

3. System overview

3.01 System specifications

Operating system Windows®10

Fast boot up time <25 seconds

Cold boot up time <100 seconds

3.02 Applications

Abdominal

Obstetrical

Gynecological

Small parts

Musculoskeletal

Vascular/peripheral vascular Urological

Pediatric

Cardiac

Thoracic

Transcranial

Transvaginal

Transrectal

Interventional Guidance

2. User interface

2.01 Operator keyboard

Full alphanumeric keypad covered with washable protection film

6 TGC pods

21.5" (476.1 x 267.8 mm) 1920 x 1080 high-resolution LED
backlit

3.03 Scanning methods

Electronic convex

Electronic linear

Electronic micro convex

Electronic sector

Mechanical volume sweep

3. System overview (cont.)

3.04 Transducer Type

Convex array
Linear array
Microconvex array
Sector phased array

Breast Productivity

Urological Calculations

Renal calculations

Cardiac calculations

On board reporting package

Network storage

Remote capability: RSVP

My Trainer

Scan Assistant

Standby

QAnalysis

3.05 Operating modes

B-Mode
Coded Harmonic Imaging
M-Mode
Color M-Mode
Color Flow mode (CFM)
Power Doppler Imaging (PDI)
Directional PDI
PW Doppler with high PRF

4. System standard features

Installation wizard
Whizz
CrossXBeam™
SRI-HD (High Definition Speckle Reduction Imaging)
B-Steer
Coded Harmonic Imaging
Virtual Convex
Patient information database
Image Archive on integrated SSD
Raw Data Analysis
Real-time automatic Doppler calculations
OB Calculations
Fetal Trending
Multi-gestational calculations
Hip dysplasia calculations
Gynecological calculations
Vascular Calculations

5. System options

CW Doppler mode

Anatomical and Curved M Mode (AMM and Curved AMM)

LOGIQView

Advanced 3D (Easy 3D)

Tissue Velocity Imaging (TVI)

TVM

Auto Bladder (Dynamic image optimization, Auto measurement
and Auto annotation)

Scan Coach

Sono Biometry (HC/ AC/ HL/ FL)

Auto IMT

TUI

Thyroid productivity A package in thyroid measurement including measurement and relevant description (Includes TI-RADS® ACR)

Needle recognition

Follow up tool

Elastography

Breast Care

Li Rads

Tricefy™ Uplink

6. Peripheral options

Sony UP-D898MD B/W thermal printer
Sony UP-D898DC B/W thermal Printer
Sony UP-D25MD Color thermal printer
1-Pedal and 3-Pedal type footswitch
USB Stick
External USB HDD
DVD RW kit
USB wireless adaptor: Sales availability varies in different countries
ECG Module
HP Office 200 Printer
Bluetooth adapter
Barcode Scanner

8. Selectable alternating modes

B/M/PW/CW/CF/PDI/TVI/TVD
B + B
B + M
B + PW/CW
B + CFM/PDI
B + CFM /PDI + PW /CW
B + TVI
B + TVI + TVD
Multi-image split screen (quad screen) Live/Frozen
Independent CINE playback
Bluetooth adapter
Barcode Scanner

7. Display modes

7.01 Live and stored display format	
Widescreen	<ul style="list-style-type: none">• Full size and split screen• Both with thumbnails for still and Cine
Review image format	4x4 and thumbnails for still and Cine
Simultaneous capability	<ul style="list-style-type: none">• Dual B (B/B)• B + PW/M• Real-time triplex mode (B + CFM/ PDI + PW)• B + CFM/PDI• B + CFM + M• B + B-Flow/ B-Flow Color
Zoom	Write (HD)/readup to 67X
Colorized Image	<ul style="list-style-type: none">• Colorized B• Colorized M• Colorized CW• Colorized 4D• Colorized B-Flow• Colorized PW• Colorized 3D
Timeline display	<ul style="list-style-type: none">• Independent dual B/PW or CW display• Display Format<ul style="list-style-type: none">- Top/bottom selectable format (Size: 1/2:1/2; 1/3:2/3; 2/3:1/3)- Side/side selectable format (Size: 1/2:1/2; 1/4:3/4; TL only)
LOGIQ View	
VirtualConvex	
TUI (Tomography Ultrasound Imaging)	

9. Display annotation

9.01 General user interface			
Patient name:	First, last	(up to 64 total characters)	
Patient ID	(Up to 64 characters)		
Other ID	(Up to 64 characters)		
Age, gender and date of birth			
Hospital name			
Date format: 4 Types selectable	<ul style="list-style-type: none">• MM/DD/YYYY• YYYY/MM/DD	<ul style="list-style-type: none">• DD/MM/YYYY• YYYY-MM-DD	
Time format: 2 types selectable	<ul style="list-style-type: none">• 24 hours• 12 hours		
Gestational age from	<ul style="list-style-type: none">• LMP• EDD	<ul style="list-style-type: none">• GA• BBT	
Displayed Acoustic	<ul style="list-style-type: none">• TIs: Thermal Index Soft Tissue• Tlc: Thermal Index Cranial (Bone)• Tib: Thermal Index Bone• MI: Mechanical Index		
% of maximum power output			
Map name			
Probe orientation			
Depth scale marker			
Lateral scale marker			
Focal zone marker			
Image depth			
Zoom depth			

9. Display annotation (cont.)

9.02 B-Mode

Gain
Dynamic range
Imaging frequency
Edge enhance
Frame average
Frame rate
Gray map
SRI-HD
CrossXBeam

9.05 M-Mode

Gain
Dynamic Range (Use the Dynamic Range of B-Mode)
Time Scale
AMM

9.06 Doppler Mode

Gain
Angle Correct
Sample Volume Depth and Sample Volume Length
Wall Filter
Baseline
Spectrum Inversion
Time Scale/Sweep Speed
Scale
Doppler Frequency

9.07 Easy 3D

Utilities
Texture
Gray Surface
Render
Threshold1
Threshold2
Scan Distance
Colorize

9.08 Advanced 3D

DefineAxis
Group Planes
Reslice
Tile

9.04 PDI mode

Line Density
Frame Average
Packet Size
Directional PDI
Color Velocity Range and Baseline
Power Threshold Marker
PDI Gain
Inversion

10. General system parameters

10.01 System Setup

10 Pre-programmable Categories

User Programmable Preset Capability

Factory Default Preset Data

Languages	English, Latin American Spanish, French, German, Italian, Brazilian Portuguese, Chinese (simplified), Swedish, Russian, Norwegian, Danish, Dutch, Finnish
-----------	---

OB Report Formats

Tokyo Univ., Osaka Univ., USA, Europe, and ASUM

User defined annotations

Body patterns

Customized comment home position

10.03 B-Mode

Acoustic power output 0 – 100%, 2, 5 and 10 steps

Gain From 0 – 90 dB, 1 dB per step

Adjustable dynamic range 36 – 96 dB, 3 or 6 dB per step

Frame averaging 8 steps

Maximum frame rate ≥ 1449 fps

Grayscale map 6 or 8 types, probe and application dependent

B colorization 9 types

Frequency Up to 4 selectable, probe dependent

Line density 992, 5–7 steps, probe dependent

Line density zoom 5–7 steps, probe dependent

Thermal index TIC, TIS, TIB

Image reverse On/off

Focus number 8 steps

Focus width 3 types

Suppression 6 steps

Edge enhance 7 steps

Rejection Up to 9 steps, probe dependent

Steered linear ±12°, ±15°, probe dependent

Scanning size (FOV or angle, probe dependent)

SRI-HD Up to 8 levels selectable

CrossXBeam Up to 9 angles selectable, probe dependent

Depth 1 – 33 cm, 0.5, 1 or 2 cm per step, probe dependent

10.02 System scanning parameters

Digital Agile Beamformer Architecture

115,663 system processing channels

Max. Frame Rate: 1790 fps, probes and modes dependent

Displayed Imaging Depth: 1 – 33 cm

Minimum Depth of Field: 0 – 1 cm, probe dependent

Maximum Depth of Field: 0 – 33 cm, probe dependent

Transmission Focus: 1 – 8 Focal Points selectable, probe and application dependent

Quad Beamforming

Continuous Dynamic Receive Focus/Aperture

Multi-Frequency/Wideband Technology

Frequency Range: 1.7 to 13 MHz

256 Shades of Gray

266 dB systematic Dynamic Range

Adjustable Field of View (FOV): Up to 168 degree, probe dependent

Image Reverse: Right/Left

Image Rotation of 0° 90° 180° 270°

10. General system parameters (cont.)

10.04 Coded Harmonic Imaging		10.07 Color Flow mode	
Available on all probes			
Line density	5 or 6 steps, probe dependent	Baseline	0 –100%, 10% per step
Line density zoom	5 or 6 steps, probe dependent	Invert	Off/on
Suppression	6 steps	CF/PDI focus depth	Default pre-settable for 10 –100% of ROI in depth, 15% or 20% per step
Edge enhance	7 steps	CF PDI flash suppression	5 steps
Gray map	6 or 8 types, probe and application dependent	CF/PDI angle steer	0, ±10°, ±15°, ±20°, probe dependent
Tint map	9 types	Packet size	8 –24, probe and application dependent
Gain	0 – 90 dB, 1 dB per step	Line density	5 steps
Dynamic range	51 –78 dB, 3 dB per step; 36 –48 dB/78 –96 dB, 6 dB per step;	Line density zoom	5 steps
Rejection	Up to 9 steps, probe dependent	Frame average	7 steps
Frequency	Up to 4 steps, probe dependent	Maximum frame rate	≥1449 fps
10.05 SRI-HD		PRF	0.1 –22.3 KHz
High Definition Speckle Reduction Imaging provides multiple levels of speckle reduction		Spatial filter	6 steps
Compatible with side-by-side DualView display		Gain	0 – 40 dB, 0.5 dB per step
Compatible with all linear, convex and sector transducers		Wall filter	4 steps, probe and application dependent
Compatible with B-Mode, 3D/4D imaging		Scanning size (FOV or angle)	Probe dependent
10.06 CrossXBeam		CF/PDI vertical size (mm) of ROI	Default pre-settable
Provides 3, 5, 7, 9 of spatial compounding		CF/PDI vertical size (mm) of ROI	Default pre-settable
Live side-by-side DualView display		CF/PDI frequency	Up to 4 steps, probe dependent
Compatible with	<ul style="list-style-type: none"> • Color Mode • PW • SRI-HD • Coded Harmonic Imaging • Virtual Convex 	Colormaps, including velocity-variance maps	19 types, probe and application dependent
Available on 4C-RS, L6-12-RS, E8C-RS, E8Cs-RS, 8C-RS, RAB2-6-RS, LK760-RS		Transparent map	5 steps
		Colorthreshold	0 –100%, 10% per step
		Accumulation	8 steps

10. General system parameters (cont.)

10.08 Power Doppler Imaging mode		10.10 Anatomical M-Mode (option)
DI map	14 types	M-Mode cursor adjustable at any plane
CF/PDI focus depth	Default pre-settable for 10 –100% of ROI in depth, 15% or 20% per step	Can be activated from a Cine loop from a live or stored image
PDI acoustic output	0 – 100%, 2%, 5% or 10% per step	Measure and analysis capability
CF/PDI angle steer	0, ±10°, ±15°, ±20°, probe dependent	Available with ColorFlow mode
Packet size	8 –24, probe and application dependent	
Spatial filter	6 steps	
Frame average	7 steps	
PRF	0.1–22.3 KHz	
Power threshold	0 –100%, 10% per step	
Gain	0 – 40 dB, 0.5 dB per step	
Wall filter	4 steps, probe and application dependent	
CF/PDI frequency	Up to 4 steps, probe dependent	
Transparent map	5 steps	
Invert	On/off	
Accumulation	8 steps	
10.09 M-Mode		10.11 Pulse Wave Doppler mode
Gain	-20 – 20 dB, 1 dB per step	Acoustic power 0 – 100%, 2, 5 and 10 steps
Grayscale map	6 or 8 types, probe dependent	Gain 0 – 85 dB, 1 dB per step
Colorization	9 types	Grayscale map Up to 8 types
Scanning size (FOV or angle, probe dependent, see probe specifications)		PRF 0.6 –27.9 KHz
Rejection	6 steps	Transmit frequency 1.7 –6.3 MHz, probe dependent
Compression	13 steps	Wall filter 5.5 –5000 Hz, 27 steps, probe dependent
Sweep Speed	8 steps	PW colorization Up to 6 types
M/PW display format	Vert 1/3B, Vert 1/2B, Vert 2/3B, Horiz 1/2B, Horiz 1/4B, TL only	Velocity scale range 0.1 –7011 cm/s
		Sample volume depth 0.1 –33 cm, probe dependent
		Sweep speed 0 –7, 8 steps
		SV gate 1, 2, 3, 4, 5, 6, 7, 8, 10, 12, 14, 16 mm
		Angle correction -90° –90°, 1° per step
		M/PW display format Vert 1/3B, Vert 1/2B, Vert 2/3B, Horiz 1/2B, Horiz 1/4B, TL only
		Spectrum inversion Off/on
		Simultaneous Off (PW only)/on
		PW angle steer 0, ±10°, ±15°, ±20° (use angle steer of B Mode), probe dependent
		Tracemethod Off, Max, Mean
		Baseline shift 11 steps
		Auto Calcs/Doppler Off, Frozen, Live
		AutoTrace
		Compression 0.5 –2.4 (0.5, 0.7, 0.9, 1, 1.1, 1.4, 1.6, 2, 2.4)
		Tracedirection Above, Below, Both
		Tracesensitivity 0 –40, 2 per step

10. General system parameters (cont.)

10.12 Continues Wave Doppler mode		10.14 Image storage
Grayscale map	8 types	On-board database of patient information
Baseline	11 steps	Conversion to formats JPEG, AVI, WMV
Angle correct	-90° –90°, 1° per step	Live image and stored image side-by-side display
Spectralcolor	6 types	Reload of archived data sets
Invert	Off/on	Network storage support for Import, Export, DICOM Read, SaveAs
Cycles to average/ Spectral averaging	5 steps	Storage formats <ul style="list-style-type: none"> DICOM – compressed/uncompressed, single/multi-frame, with/without Raw Data Export JPEG, WMV (MPEG 4) and AVI formats DICOM still image storage size: ~3.9 MB Display format: Full size, 4x4 and thumbnails
Gain	0 – 85 dB, 1 dB per step	Storage devices <ul style="list-style-type: none"> Internal Solid State drive partition of 100 GB for image storage Optional internal SSD of 256GB for image storage External USB HDD and USB memory stick support for Import, Export, DICOM Read, SaveAs CD-RW storage: 700 MB DVD storage: -R (4.4 GB)
Wall filter	5.5 –5000 Hz, 27 steps, probe and application dependent	
CW-Mode includes	<ul style="list-style-type: none"> Transmit frequency: 1.9, 4.2 MHz CW colorization: Tint map A/B/C/D/E/F Velocity scale range: 0.2 – 6105 cm/s Spectrum inversion Trace method: Max, Mean, Off Auto Calcs/Doppler Auto Trace: Frozen, Live, Off Trace direction: Above, Below, Both 	
Trace sensitivity:	0 – 40, 2 per step	
10.13 Cine memory/image memory		10.15 Connectivity and DICOM
384 MB of Cine memory		Ethernet network connection <ul style="list-style-type: none"> DICOM 3.0 Verify Print Store Modality worklist Storage commitment Modality Performed Procedure Step (MPPS) Query/retrieve Structured reporting template (Can be compared to vascular and OB standard) RSVP Remote capability
Selectable Cine sequence for Cine review		
Maximum number of cine loops	410691 images	
Cine loop capacity	<ul style="list-style-type: none"> B-Mode: 919 seconds M-Mode: 423 seconds CFM Mode: 4978 seconds PW Mode: 27983 seconds CW Mode: 27962 seconds 	
Prospective Cine mark		
Measurements/calculations and annotations on Cine playback		
Scrolling timeline memory		
Dual image Cine display		
Quad image Cine display		
Cine gauge and Cine image number display		
Cine review loop		
Cine review speed	11 steps (11, 13, 14, 17, 22, 25, 31, 48, 100, 200, 400%)	Provides a convex field of view
Compatible with CrossXBeam for linear transducers		
Available on linear and sector transducers		

10. General system parameters (cont.)

10.17 LOGIQView (Option)

Extended Field of View Imaging

Available on 4C-RS, L6-12-RS, 8C-RS, 3Sc-RS, E8C-RS, E8Cs-RS, 6S-RS, LK760-RS, RAB2-6-RS probes

Foruse in B-Mode

CrossXBeam is available on linear probes

Auto detection of scan direction

Post-processz oom Rotation

Auto fit on monitor

Measurements in B-Mode

Up to 60 cm scan length

10.18 Easy 3D (Option)

Allows unlimited rotation and planar translations

3D reconstruction from Cine sweep

Threshold1: 0 –255

Utilities: Average off/Average light/Average medium/Average Strong

Grey surface: Off/On

Scan distance: 1.0 –15.0

Threshold2: 0 –255

Colorize: 0 –360

10.19 Advanced 3D (Option)

Define axis: Select 2 points as start and end point of long axis

Group planes: Off/Main/Parallel/Angular

Reslice: Cube/Virtual Rescan/Cubic Plane

Tile: 1/2/4/6

10.20 TVI (option)

Myocardial Doppler imaging with color overlay on tissue image

Available on the sector probes

Tissue color overlay can be removedto show just the 2D image, still retaining the tissue velocity information

QAnalysis: Multiple Time Motion trace display from selected points in the myocardium

10.21 TVM (option)

TVI with M-Mode active

Available on the sector probes

Provides both myocardium motion velocity and direction

10.22 Follow-up tool (option)

The follow-up tool is intended to more accurately perform serial scans on a patient, and compare the images of a previous ultrasound exam with the current exam

10.23 Needle recognition (option)

Needle recognition allows you to obtain precise needle imaging in the dashed box. It is available with probes on L6-12-RS, 4C-RS.

10.24 Scan Coach (option)

Scan Coach is a contextual reference tool. It is with clinical guidance for scan plane acquisition and references for anatomical structures. It can be displayed on-demand by the user. Clinical reference images and animations to depict information related to each step. It covers five applications.

- Abdomen
- Obstetrics
- Gynecology
- Cardiology
- Vascular

10.25 My Trainer

Abstracted from basic user manual, it lists out FAQs from customers and instructs customer how to solve problems by themselves timely.

10.26 Battery (option)

The lithium ion battery provides power when an AC power source is not available. About 15 minutes of battery life can be expected with fully charged battery in use to supply power to the system.

10.27 Scan Assistant

Scan Assistant provides an automated exam script that moves you through an exam step-by-step. This allows you to focus on performing the exam rather than on controlling the system and can help you to increase consistency while reducing keystrokes.

10. General system parameters (cont.)

10.28 RSVP

RSVP is a direct link with a GE Online Service Engineer or Applications Support Engineer or a Request for Service.

10.29 Whizz

Whizz will continuously optimize the brightness, contrast and uniformity of B-Mode images when scanning different tissues. Whizz in PW/CW Doppler Mode optimizes the spectral data. Auto adjusts the Velocity Scale/PRF (live imaging only), baseline shift, and invert (if preset). Upon deactivation, the spectrum is still optimized.

10.30 Whizz CF mode

Whizz CF mode dynamically optimizes CFmode parameters in real time, enabling optimal, consistent image quality and measurements.

10.31 Lateral gain compensation (LGC)

To set Lateral Gain Compensation values based on LGC curves defined by a user; By laterally adjusting the received signal intensity, the uniformity of a B-Mode image intensity is optimized in lateral direction.

10.32 Elastography (option)

Available on L6 12 RS

- Frame reject: 0 – 8
- Axial smoothing: 0 – 4
- Noise reject: 0 – 8
- Sample Volume: 0 – 2
- Lateral smoothing: 0 – 4
- Window: 0 – 8
- Map: 8
- Frame average: 0 – 10
- Line Density: 0 – 4
- Soft compress: 0 – 10
- Hard compress: 0 – 10

10.33 LI-RADS® ACR (Option)

The Ultrasound Liver Imaging Reporting And Data System (US LI-RADS) is a standardized system for imaging technique, interpretation, reporting, and data collection for screening or surveillance ultrasound exams in patients at risk for developing hepatocellular carcinoma(HCC).

10.34 e-Delivery

Electronic software delivery. As part of the product lifecycle management, GE regularly analyzes and integrates software updates into our products.

10.35 Probe Check

Probe assessment tool that evaluates probe elements to monitor potential probe deterioration over probe life cycle.

10.36 Controls available while “live”

Write Zoom

B/M-Mode	<ul style="list-style-type: none">• Gain• Dynamic Range• Transmission Focus Position• Line Density Control• Sweep Speed for M-Mode	<ul style="list-style-type: none">• TGC• Acoustic Output• Transmission Focus Number• Number of Angles for CrossXBeam
PW-Mode	<ul style="list-style-type: none">• Gain• Acoustic Output• Transmission Frequency• Scale• Wall Filter	<ul style="list-style-type: none">• Sample Volume Gate– Length– Depth• Volume
Color Flow Mode	<ul style="list-style-type: none">• CFM Gain• Acoustic Output• Wall Filter• Line Density• CFM Frame Average	<ul style="list-style-type: none">• CFM Velocity Range• Packet Size• CFM Spatial Filter• Frequency/Velocity Baseline Shift

10.37 Controls available on Freeze or Recall

SRI-HD

CrossXBeam –Display non-compounded and compounded image simultaneously in split screen

Easy 3D reconstruction from a stored Cine loop

CrossXBeam is disabled on Freeze or Recall TGC

Colorized B and M

Frame average (loops only)

Dynamic range

Anatomical M-Mode

Gray Map

Post gain

10. General system parameters (cont.)

10.37 Controls available on Freeze or Recall (cont.)

Baseline shift (PW, CW)

Sweep speed

Compression

Rejection

Colorized spectrum

Display format

Angle Correct

Quick Angle Correct

Overall gain (loops and stills)

Color map

Transparency map

CFM display threshold

Invert for Color/Doppler

11.03 General Doppler measurements/calculations

Velocity

Time

A/B ratio(velocities)

PS (Peak Systole)

ED (End Diastole)

PS/ED (PS/ED ratio)

ED/PS (ED/PS ratio)

AT (Acceleration Time)

ACCEL(Acceleration)

TAMAX (Time Averaged Maximum Velocity)

Volume Flow (TAMEAN and vessel area)

Heart Rate

PI (Pulsatility Index)

RI (Resistivity Index)

11. Measurements/calculations

11.01 General B-Mode

Depth and distance

Circumference (ellipse/trace)

Area (ellipse/trace)

Volume

% Stenosis (areaor diameter)

Angle between 2 lines

11.04 Real-time Doppler Auto measurements/calculations

PS (Peak Systole)

ED (End Diastole)

MD (Minimum Diastole)

PI (Pulsatility Index)

RI (Resistivity Index)

AT (Acceleration Time)

ACC (Acceleration)

PS/ED (PS/ED Ratio)

ED/PS (ED/PSRatio)

HR (Heart Rate)

TAMAX (Time Averaged Maximum Velocity)

PVAL (Peak Velocity Value)

Volume Flow (TAMEAN and Vessel Area)

11.02 General M-Mode

M-Depth

Distance

Time

Slope

Heart rate

11. Measurements/calculations

11.05 OB measurements/calculations

- | | |
|---------------------------------|--|
| Gestational age by | <ul style="list-style-type: none">• GS (Gestational Sac)• CRL (Crown Rump Length)• FL (Femur Length)• BPD (Biparietal Diameter)• AC (Abdominal Circumference)• HC (Head Circumference)• APTD x TTD (Anterior/Posterior Trunk Diameter by Transverse Trunk Diameter) |
| Trunk Diameter | <ul style="list-style-type: none">• FTA (Fetal Trunk Cross-sectional Area)• HL (Humerus Length)• BD (Binocular Distance)• FT Foot Length)• OFD (Occipital Frontal Diameter)• TAD Transverse Abdominal Diameter)• TCD (Transverse Cerebellum Diameter)• THD (Thorax Transverse Diameter)• TIB (Tibia Length)• ULNA (Ulna Length) |
| Estimated Fetal Weight (EFW) by | <ul style="list-style-type: none">• AC, BPD• AC, FL, HC |
| Calculations and ratios | <ul style="list-style-type: none">• FL/BPD• FL/HC• CI (Cephalic Index) |
| SonoBiometry | <ul style="list-style-type: none">• BPD• AC• FL |

Measurements/calculations by: ASUM, ASUM 01, Berkowitz, Bertagnoli, Brenner, Campbell, CFEF, Chitty, Eik-Nes, Eriksen, Goldstein, Hadlock, Hansmann, Hellman, Hill, Hohler, Jeanty, JSUM, Kurtz, Mayden, Mercer, Merz, Moore, Nelson, Osaka, Paris, Rempen, Robinson, Shepard, Shepard/Warsoff, Tokyo, Tokyo/Shinozuka, Yarkoni

Fetal graphical trending

Growth percentiles

Multi gestational calculations

Fetal qualitative description (anatomical survey)

Fetal Environmental Description (Biophysical profile)

Programmable OB tables

Over 20 selectable OB calculations

Expanded worksheets

11.06 GYN measurements/calculations

Right ovary length, width, height

Left ovary length, width, height

Uterus length, width, height

Cervix length, trace

Ovarian volume

ENDO (Endometrial thickness)

Ovarian RI

Uterine RI

Follicular measurements

11.07 Vascular measurements/calculations

DCCA (Distal Common Carotid Artery)

MCCA (Mid Common Carotid Artery)

PCCA (Proximal Common Carotid Artery)

DICA (Distal Internal Carotid Artery)

MICA (Mid Internal Carotid Artery)

PICA (Proximal Internal Carotid Artery)

DECA (Distal External Carotid Artery)

PECA (Proximal External Carotid Artery)

VERT (Vertebral Velocity)

SUBCLAV (Systolic Subclavian Velocity)

Automatic IMT

11.08 Urological calculations

Volume (Auto Bladder volume) Prostate volume

Left/right renalvolume

Generic volume

Post-void bladder volume

12. Cardiac measurements/calculations

12.01 B-Mode measurements		12.02 M-Mode measurements	
Aorta	<ul style="list-style-type: none"> • Aortic Root Diameter (Ao Root Diam) • Aortic Arch Diameter (Ao Arch Diam) • Ascending Aortic Diameter (Ao Asc) • Descending Aortic Diameter (Ao Desc Diam) • Aorta Isthmus (Ao Isthmus) • Aorta (Aost junct) 	Aorta	<ul style="list-style-type: none"> • Aortic Root Diameter (Ao Root Diam) • Aortic Valve Diameter AV Diam) • Aortic Valve Cusp Separation AV Cusp) • Aortic Valve Ejection Time (LVET)
Aortic valve	<ul style="list-style-type: none"> • Aortic Valve Cusp Separation (AV Cusp) • Aortic Valve Area Planimetry (AVA Planimetry) • (Trans AVA) 	Left atrium & Left ventricle	<ul style="list-style-type: none"> • Left Atrium Diameter to AoRoot Diameter Ratio (LA/Ao Ratio) • Left Atrium Diameter (LA Diam) • Left Ventricle Volume, Teichholz/Cubic (LVIDd, LVI Ds) • Left Ventricle Posterior Wall Thickness (LVPWs) • Left Ventricle Ejection Time (LVET) • Left Ventricle Pre Ejection Period (LVPEP) • Interventricular Septum (IVS)
Left atrium	<ul style="list-style-type: none"> • Left Atrium Diameter (LA Diam) • LA Length (LA Major) • LA Width (LA Minor) • Left Atrium Area (LAA(d), LAA(s)) • Left Atrium Volume, Single Plane, Method of Disk (LAEDV A2C, LAESV A2C) (LAEDV A4C, LAESV A4C) 	Mitral valve	<ul style="list-style-type: none"> • E-Point-to-Septum Separation (EPSS) • Mitral Valve Anterior Leaflet Excursion (D-E Excursion) • Mitral Valve D-E Slope (D-E Slope) • Mitral Valve E-F Slope (E-F Slope)
Left ventricle	<ul style="list-style-type: none"> • Left Ventricle Volume, Teichholz/Cubic (LVIDd, LVIDs) • Left Ventricle Internal Diameter (LVIDd, LVIDs) • Left Ventricle Length (LVLD, LVLs) • Left Ventricle Outflow Tract Diameter (LVOT Diam) • Left Ventricle Posterior Wall Thickness (LVPWd, LVPWs) • Left Ventricle Length (LV Major) • Left Ventricle Width (LV Minor) • Left Ventricle Outflow Tract Area (LVOT) • Left Ventricle Mass Index (LVPWd, LVPWs) • Ejection Fraction, Teichholz/Cube (LVIDd, LVIDs) • Left Ventricle Posterior Wall Fractional Shortening (LVPWd, LVPWs) • Mitral Valve • E-Point-to-Septum Separation (EPSS) • Mitral Valve Area Planimetry (MVA Planimetry) 	Pulmonic valve	<ul style="list-style-type: none"> • QRS complex to end of envelope (Q-to-PV close) • Right Ventricle Internal Diameter (RVIDd, RVIDs) • Right Ventricle Outflow Tract Diameter (RVOT Diam)) • Right Ventricle Ejection Time (RVET) • Right Ventricle Pre-Ejection Period (RVPEP)
Pulmonic valve	Pulmonic Diameter (Pulmonic Diam)	Tricuspid valve	<ul style="list-style-type: none"> • QRS complex to end of envelope (Q-to-TV close)
Right ventricle	<ul style="list-style-type: none"> • Right Ventricle Internal Diameter (RVIDd, RVIDs) • Right Ventricle Outflow Tract Diameter (RVOT Diam) 	Aortic valve	<ul style="list-style-type: none"> • Aortic Valve Mean Velocity (AV Trace) • Aortic Valve Velocity Time Integral (AV Trace) • Aortic Valve Mean Pressure Gradient (AV Trace) • Aortic Valve Peak Pressure Gradient (AR Vmax) • Aortic Insufficiency Peak Velocity (AR Vmax) • Aortic Insufficiency End-Diastolic Velocity (AR Trace) • Aortic Valve Peak Velocity (AV Vmax) • Aortic Valve Deceleration Time (AV Trace) • Aortic Valve Ejection Time (AVET) • Aortic Valve Area according to PHT
System inferior vena cava	<ul style="list-style-type: none"> • Systemic Vein Diameter (Sytemic Diam) 		

12.03 Doppler mode measurements

Aortic valve	<ul style="list-style-type: none"> • Aortic Valve Mean Velocity (AV Trace) • Aortic Valve Velocity Time Integral (AV Trace) • Aortic Valve Mean Pressure Gradient (AV Trace) • Aortic Valve Peak Pressure Gradient (AR Vmax) • Aortic Insufficiency Peak Velocity (AR Vmax) • Aortic Insufficiency End-Diastolic Velocity (AR Trace) • Aortic Valve Peak Velocity (AV Vmax) • Aortic Valve Deceleration Time (AV Trace) • Aortic Valve Ejection Time (AVET) • Aortic Valve Area according to PHT
--------------	--

12. Cardiac measurements/calculations (cont.)

12.03 Doppler mode measurements (cont.)		12.03 Doppler mode measurements (cont.)
Left ventricle	<ul style="list-style-type: none"> • Left Ventricle Outflow Tract Peak Pressure Gradient (VLOT Vmax) • Left Ventricle Outflow Tract Peak Velocity (LVOT Vmax) • Left Ventricle Outflow Tract Mean Pressure Gradient (LVOT Trace) • Left Ventricle Outflow Tract Velocity Time Integral (LVOT Trace) • Left Ventricle Ejection Time (LVET) 	<ul style="list-style-type: none"> • Pulmonic Valve Peak Velocity (PV Vmax) • Pulmonary Artery Diastolic Pressure (PV Trace) • Pulmonic Insufficiency Mean Pressure Gradient (PR Trace) • Pulmonic Valve Mean Pressure Gradient (PV Trace) • Pulmonic Insufficiency Mean Square Root Velocity (PR Trace) • Pulmonic Insufficiency Velocity Time Integral (PR Trace) • Pulmonic Valve Mean Velocity (PV Trace) • Pulmonic Valve Velocity Time Integral (PV Trace) • Pulmonic Insufficiency Pressure Half Time (PR PHT) • Pulmonic Valve Flow Acceleration (PV Acc Time) • Pulmonic Valve Acceleration Time (PV Acc Time) • Pulmonic Valve Ejection Time (PVET) • QRS complex to end of envelope (Q-to-PV close) • Pulmonic Valve Acceleration to Ejection Time Ratio (PV Acc Time, PVET)
Mitral valve	<ul style="list-style-type: none"> • Mitral Valve Regurgitant Mean Velocity (MR Trace) • Mitral Regurgitant Mean Pressure Gradient (MR Trace) • Mitral Regurgitant Velocity Time Integral (MR Trace) • Mitral Valve Mean Velocity (MR Trace) • Mitral Valve Velocity Time Integral (MR Trace) • Mitral Valve Mean Pressure Gradient (MR Trace) • Mitral Regurgitant Peak Pressure Gradient (MR Vmax) • Mitral Valve Peak Pressure Gradient (MR Vmax) • Mitral Regurgitant Peak Velocity (MR Vmax) • Mitral Valve Peak Velocity (MR Vmax) • Mitral Valve Velocity Peak A (MV A Velocity) • Mitral Valve Velocity Peak E (MV E Velocity) • Mitral Valve Area according to PHT (MV PHT) • Mitral Valve E-Peak to A-Peak Ratio (A-C and D-E) (MV E/ARatio) • Mitral Valve Acceleration Time (MV ACC Time) • Mitral Valve Deceleration Time (MV Dec. Time) • Mitral Valve Acceleration Time/Deceleration Time Ratio (MVAcc/Dec. Time) 	<ul style="list-style-type: none"> • Right Ventricle Outflow Tract Peak Pressure Gradient (RVOT Vmax) • Right Ventricle Outflow Tract Peak Velocity (RVOT Vmax) • Right Ventricle Outflow Tract Velocity Time Integral (RVOT Trace) • Right Ventricle Ejection Time (RV Trace) • Stroke Volume by Pulmonic Flow (RVOT Planimetry, RVOT Trace) • Right Ventricle Stroke Volume Index by Pulmonic Flow (RVOT Planimetry, RVOT Trace)
Pulmonic valve	<ul style="list-style-type: none"> • Pulmonic Insufficiency Peak Pressure Gradient (PR Vmax) • Pulmonic Insufficiency End-Diastolic Pressure Gradient (PRTrace) • Pulmonic Valve Peak Pressure Gradient (PV Vmax) • Pulmonic Insufficiency Peak Velocity (PR Vmax) • Pulmonic Insufficiency End-Diastolic Velocity (Prend Vmax) 	<ul style="list-style-type: none"> • Pulmonary Artery Peak Velocity (PV Vmax) • Pulmonary Vein Velocity Peak A (reverse) (P Vein A) • Pulmonary Vein Peak Velocity (P Vein D, P Vein S) • Systemic Vein Peak Velocity (PDA Diastolic, PDA Systolic) • Ventricular Septal Defect Peak Velocity (VSD Vmax) • Atrial Septal Defect (ASD Diastolic, ASD Systolic) • Pulmonary Vein A-Wave Duration (P Vein A Dur) • IsoVolumetric Relaxation Time (IVRT) • IsoVolumetric Contraction Time (IVCT) • Pulmonary Vein S/D Ratio (P Vein D, P Vein S) • Ventricular Septal Defect Peak Pressure Gradient (VSD Vmax) • Pulmonic-to-Systemic Flow Ratio (Qp/Qs)

12. Cardiac measurements/calculations (cont.)

12.03 Doppler mode measurements (cont.)		12.05 Combination mode measurements	
Tricuspid valve	<ul style="list-style-type: none"> • Tricuspid Regurgitant Peak Pressure Gradient (TR Vmax) • Tricuspid Valve Peak Pressure Gradient (TVVmax) • Tricuspid Regurgitant Peak Velocity (TR Vmax) • Tricuspid Valve Peak Velocity(TVVmax) • Tricuspid Valve Velocity Peak A (TV A Velocity) • Tricuspid Valve Velocity Peak E (TV E Velocity) • Tricuspid Regurgitant Mean Pressure Gradient (TRTrace) • Tricuspid Valve Mean Pressure Gradient (TVTrace) • Tricuspid Regurgitant Velocity Time Integral (TR Trace) • Tricuspid Valve Mean Velocity (TV Trace) • Tricuspid Valve Velocity Time Integral (TV Trace) • Tricuspid Valve Time to Peak (TV Acc/Dec Time) • Tricuspid Valve Ejection Time (TV Acc/Dec Time) • Tricuspid Valve A-Wave Duration (TVA Dur) • QRS complex to end of envelope (Q-to-TV close) • Tricuspid Valve Pressure Half Time (TV PHT) • Tricuspid Valve E-Peak to A-Peak Ratio (TV E/A Velocity) 	Aortic valve	<ul style="list-style-type: none"> • Aortic Valve Area (Ao Diam., LVOT Vmax, AV Vmax) • Aortic Valve Area by Continuity Equation by Peak Velocity (Ao Diam, LVOT Vmax, AV Vmax) • Stroke Volume by Aortic Flow (AVA Planimetry, AV Trace) • Cardiac Output by Aortic Flow (AVA Planimetry, AV Trace, HR) • Aortic Valve Area by Continuity Equation VTI (Ao Diam, LVOT Vmax, AV Trace)
		Left ventricle	<ul style="list-style-type: none"> • Cardiac Output, Teichholz/Cubic (LVIDd, LVI Ds, HR)
		Mitral valve	<ul style="list-style-type: none"> • Stroke Volume by Mitral Flow (MVA Planimetry, MV Trace) • Cardiac Output by Mitral Flow (MVA Planimetry, MV Trace, HR)
12.04 Color Flow mode measurements		12.06 Cardiac worksheet	
Aortic valve	<ul style="list-style-type: none"> • Proximal Isovelocity Surface Area: Regurgitant Flow (AR Trace) • Proximal Isovelocity Surface Area: Regurgitant Volume Flow (AR Trace) • Proximal Isovelocity Surface Area: Aliased Velocity (AR Vmax) 	<p>Parameter: Lists the mode, the measurement folder and the specific measurement</p> <p>Measured Value: Up to six measurement values for each item. Average, maximum, minimum, or last</p> <p>Generic study in cardiology</p>	
Mitral valve	<ul style="list-style-type: none"> • Proximal Isovelocity Surface Area: Regurgitant Flow (MR Trace) • Proximal Isovelocity Surface Area: Regurgitant Volume Flow (MR Trace) • Proximal Isovelocity Surface Area: Aliased Velocity (MR Vmax) 		

13. Probes

13.01 4C-RS

Convex probe

Applications	Abdominal, OB, GYN, Vascular, Urology, Thoracic, Pediatric, MSK, Interventional Guidance
Number of elements	128
Convex radius	60 mm
FOV	58°
Footprint	66.2 x 18.3 mm
B-Mode imaging frequency	2.0, 3.0, 4.0, 5.0 MHz
Harmonic imaging frequency	3.0, 4.0, 5.0 MHz
CFM/PDI/PWD frequency	2.0 MHz (CFM/PDI) 2.5, 2.8, 3.3 MHz
Biopsy guide	Multi-angle, reusable bracket

13.03 LK760-RS

Linear probe

Applications	MSK, Interventional Guidance
Number of elements	128
Footprint	67.0 x 13.0 mm
B-Mode imaging frequency	5.0, 7.0, 9.0 MHz
Harmonic imaging frequency	6.0, 8.0, 10.0 MHz
CFM/PDI/PWD frequency	3.5, 4.2, 5.0 MHz
Steered angle	±12°
Biopsy guide	Not available
PWD frequency	4.0, 4.5, 5.0 MHz
Steered angle	±20°
Biopsy guide	Multi-angle, reusable bracket

13.02 L6-12-RS

Linear probe

Applications	Vascular, Pediatric, Small Parts, MSK, Thoracic, Interventional Guidance
Number of elements	128
Footprint	47 x 11.4 mm
Convex radius	60 mm
FOV	58°
B-Mode imaging frequency	6.0, 8.0, 10.0, 11.0 MHz
Harmonic imaging frequency	8.0, 10.0, 12.0, 13.0 MHz
CFM/PDI frequency	4.0, 5.0, 6.0 MHz
PWD frequency	4.0, 4.5, 5.0 MHz
Steered angle	±20°
Biopsy guide	Multi-angle, reusable bracket

13.04 E8C-RS

Endo micro convex probe

Applications	OB, GYN, Urology, Transvaginal, Transrectal, Interventional Guidance
Number of elements	128
Convex radius	10.73 mm
FOV	128°
Footprint	16.9 x 21.2 mm
B-Mode imaging frequency	6.0, 8.0, 10.0 MHz
Harmonic imaging frequency	7.0, 8.0, 10.0 MHz
CFM/PDI/PWD frequency	4.2, 5.0, 6.3 MHz
Biopsy guide	Fixed angle, disposable or reusable bracket

13. Probes (cont.)

13.05 8C-RS

Micro convex probe

Applications	Pediatric, MSK, Cardiac Pediatric, Transcranial, Interventional Guidance
Number of elements	128
Convex radius	10.73 mm
FOV	131°
Footprint	22.0 x 12.0 mm
B-Mode imaging frequency	6.0, 8.0, 10.0 MHz
Harmonic imaging frequency	6.0, 7.0, 8.0, 10.0 MHz
CFM/PDI/PWD frequency	4.2, 5.0, 6.3 MHz
Biopsy guide	Not available

13.07 6S-RS

Phased array sector probe

Applications	Cardiac Pediatric, Vascular, Pediatric, Transcranial, Interventional Guidance
Number of elements	64
FOV	120°
Footprint	23.5 x 16.8 mm
B-Mode imaging frequency	4.0, 5.0, 6.0 MHz
Harmonic imaging frequency	4.0, 5.0, 6.0, 7.0 MHz
CFM/PDI/PWD frequency	2.5 (CFM/PDI), 3.0, 4.0, 4.5 MHz
CWD frequency	4.2 MHz
Biopsy guide	Not available

13.06 3Sc-RS

Phased array sector probe

Applications	Cardiac, Abdominal, Vascular, Transcranial, Thoracic, Interventional Guidance
Number of elements	64
Footprint	23.7 x 18.4 mm
B-Mode imaging frequency	2.0, 3.0, 4.0 MHz
Harmonic imaging frequency	3.0, 3.2, 3.5, 4.0 MHz
CFM/PDI/PWD frequency	1.7, 2.0, 2.5, 3.3 MHz
CWD frequency	1.9 MHz
Biopsy guide	Multi-angle, reusable bracket
Biopsy guide	Not available

13.09 E8Cs-RS

Endo micro convex probe

Applications	OB, GYN, Urology, Transvaginal, Transrectal, Interventional Guidance
Number of elements	128
Convex radius	8.73 mm
Footprint	18.6 x 13.9 mm
FOV	168°
B-Mode imaging frequency	6.0, 8.0, 10.0 MHz
Harmonic imaging frequency	7.0, 8.0, 10.0 MHz
CFM/PDI/PWD frequency	4.0, 5.0, 6.0 MHz
Biopsy guide	Fixed angle, disposable or reusable bracket

14. Inputs and outputs

HDMI output (1920 x 1080 resolution)

S-Video output (via optional video adapter)

CVBS output (via optional video adapter) Ethernet (RJ45)

USB-A (2x in rear, 1 beside keyboard)

USB-C (1x in rear)

15. Safety conformance

The Versana Essential is CE marked to Council Directive 2017/745 on medical devices

Conforms to the following standards for safety

- IEC 60601-1 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
- IEC 60601-1-2 Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – requirements and tests EMC Emissions Group 1 Class A device requirements as per CISPR 11
- IEC 60601-2-37 Medical electrical equipment – Part 2-37: Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment
- ISO 10993-1 Biological evaluation of medical devices – Part 1 Evaluation and testing within a risk management process
- EN 62366-1 Medical devices –Part 1: Application of usability engineering to medical devices

