







Versana Premier™ Lotus

Transducer guide





Featuring VisionBoost Technology

Versana Premier Lotus is designed to help you apply your clinical and imaging skills to deliver high-quality care to the wide range of patients you see every day. Its versatility covers a broad range of examinations, including abdominal, OB/GYN, cardiac, small parts, urology, pediatrics, MSK, thoracic, and more.



	Applications	Footprint	Frequency Range	FOV	Biopsy Guide
Convex					
 C1-5-RS	Abdomen, OB/GYN, Thoracic/ Pleural, MSK, Vascular, Urology, Pediatrics, Veterinary	17.2 x 69.3 mm	2.0 – 5.0 MHz	70°	H40432LE
 4C-RS	Abdomen, OB/GYN, Thoracic/ Pleural, MSK, Vascular, Urology, Pediatrics, Veterinary	18.3 x 66.2 mm	2.0 – 5.0 MHz	58°	E8385NA
 RAB2-6-RS	Abdomen, OB/GYN, Urology, Veterinary	62.2 x 34.0 mm	2.0 – 6.0 MHz	70° (B) 70° x 84° (Volume scan)	H48681ML
 8C-RS	Cardiac, MSK, Transcranial, Pediatrics, Veterinary	12.0 x 22.0 mm	4.2 – 10.0 MHz	131°	Not available
Linear					
 L4-20t-RS	Thoracic/Pleural, Small parts, MSK Conventional, MSK Superficial, Interventional Guidance, Vascular/ Peripheral Vascular, Pediatric	48.43 x 11.56 mm	5.0 – 20.0 MHz	38.4 mm	H45201BLF
 L8-18i-RS	Small Parts, MSK, Vascular, Veterinary	34.8 x 11.1 mm	6.7 – 18 MHz	25.2 mm	Not available
 L3-12-RS	Thoracic/Pleural, Small Parts, MSK Conventional, MSK Superficial, Interventional Guidance, Vascular/Periperal Vascular, Pediatrics	12.9 x 60.9 mm	3.6 – 12 MHz	51 mm	H48302AA
 12L-RS	Small Parts, Thoracic/Pleural, MSK, Vascular, Pediatrics, Veterinary	12.7 x 47.1 mm	4.2 – 13.0 MHz	38.4 mm	H48392LT H48392LL H40432LC

	Applications	Footprint	Frequency Range	FOV	Biopsy Guide
Linear (cont.)					
 9L-RS	Abdomen, Thoracic/ Pleural, Small Parts, MSK Conventional, MSK Superficial, Interventional Guidance, Vascular/Peripheral Vascular, Pediatrics	14.1 x 53 mm	4.0 – 10.0 MHz	44 mm	H4906BK
 L6-12-RS	Thoracic/Pleural, Small Parts, MSK, Vascular, Pediatrics	11.4 x 47.0 mm	4.0 – 13.0 MHz	38.4 mm	H40432LC
 LK760-RS	MSK, Veterinary	67 x 13mm	3.5 – 10.0 MHz	59.9 mm	Not available
Sector					
 12S-RS	Cardiac, Vascular, Pediatrics, Transcranial, Veterinary	17.6 x 13.2 mm	4.5 – 12 MHz	120°	Not available
 6S-RS	Cardiac, Transcranial, Vascular, Pediatrics, Veterinary	23.5 x 16.8 mm	3.0 – 7.0 MHz	120°	Not available
 3sC-RS	Abdomen, Thoracic/Pleural, Cardiac, Transcranial, Vascular, Veterinary	23.7 x 18.4 mm	1.7 – 4.0 MHz	120°	H46222LC
 6Tc-RS	Cardiac Adult, Interventional Guidance, Transesophageal	Tip (L x W x H): 45 x 14 x 12 mm	3.0 – 8.0 MHz	90°	Not available

Applications	Footprint	Frequency Range	FOV	Biopsy Guide	
Endocavitary					
 RIC5-9A-RS	OB/GYN, Transvaginal, Urology, Transrectal	20.96 x 23.35mm	4.2 – 10.0 MHz	179° (B) 179° x 120° (Volume angle)	H48691Z H46721R
 IC9-RS	Fetal/OB, GYN, Interventional Guidance, Urology, Transvaginal, Transrectal	19.6 x 13.6 mm	3.6 – 9.0 MHz	179°	H48691YX H48691YW
 E8Cs-RS	OB/GYN, Transvaginal, Urology, Transrectal	18.6 x 13.9 mm	4.0 – 10.0 MHz	168°	E8385MJ H40412LN
 E8C-RS	OB/GYN, Transvaginal, Urology, Transrectal	16.9 x 21.2 mm	4.2 – 10.0 MHz	128°	E8385MJ H40412LN
 E7C8L-RS	Urology, Transrectal	Micro Curve Part 17.2 x 10.3 mm Linear Part 57.6 x 9.1 mm	5.0 – 10.0 MHz	Micro Curve Part 136° Linear Part 51.2 mm	H40202E
 BE9CS-RS	Urology, Transrectal	19 x 18.7 mm	4.0 – 10.0 MHz	133°	E8387M

Versana Premier LOTUS and Versana Premier Lotus are different marketing representations for the same product.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

© 2025 GE HealthCare.

Versana Premier is a trademark of GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

April 2025
DOC2993174



GE HealthCare